	<p align="center">Document name: GEN-2-043 Attachment N-Amniotic Fluid AChE Only Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059815 NTD Labs SOP ID: GEN-2-043 (Att N) Revision:1</p>
-----------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Reporting Title: Amniotic Fluid AChE

Test Definition: AFT0002
Testing Location: Melville, NY
Reporting Location: Melville, NY

Description:

When an ONTD is suspected based upon maternal serum alpha-fetoprotein (AFP) screening results or diagnosed via ultrasound, analysis of AFP and acetylcholinesterase (AChE) in amniotic fluid are useful diagnostic tools. AChE is primarily active in the central nervous system with large amounts of non-specific cholinesterase's found in erythrocytes, skeletal muscle, and fetal serum. Normal amniotic fluid does not contain the specific neurally derived AChE, except in the presence of a fetal ONTD.

Analytical Method(s):

Acetylcholinesterase is evaluated using Gel Eletrophoresis.

Patient preparations:

Prepare patient for amniocentesis in accordance with office procedure.

Specimen Requirements:

Container/Tube: Green Top (Sodium Heparin) Vacutainer® Tube / Transfer tube

Specimen Volume: 1 mL supernatant

Specimen Stability: Amniotic fluid supernatant samples are stable at ambient temperature for 6 days.

Specimen Rejection Criteria: Incorrect tube type

Additional Information:

1. Indications for Testing: Suspicion of Open neural tube defect due to elevated Amniotic Fluid AFP or ultrasound findings.
2. Special Timing: None

CPT Code: 1 x 82013

Reference Values:

Negative (reported as negative [normal] or positive [abnormal] for inhibitable acetylcholinesterase)

An Interpretive Report will be provided.

Supplemental Report:


No

Testing Algorithm:

Follow up:

It is the experience of this program that amniotic fluids from pregnancies affected with ONTDs produce positive AChE results . Malformations of the abdominal wall of the fetus (gastroschisis or omphalocele)

Revision: 1 © This document is copyright of Eurofins	Effective date: May 18, 2018 During the actual use of this form, cover page need not be printed.	Page 1 of 8
---------------------------------------------------------	------------------------------------------------------------------------------------------------------------	-------------

	<p align="center">Document name: GEN-2-043 Attachment N-Amniotic Fluid AChE Only Instruction Manual</p>	<p align="center">Eurofins Document Reference: 1-D-QM-CF -9059815 NTD Labs SOP ID: GEN-2-043 (Att N) Revision:1</p>
-----------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

produce positive AChE in approximately 60% of affected cases. AChE results have been negative in cases of congenital nephrosis. Patients should be counseled accordingly.

Consents/Authorizations:

Patient signature on patient authorization/assignment on requisition form is required.

Disclaimer:

The test was developed and its performance characteristics determined by Eurofins NTD, LLC. It has not been cleared or approved by the U.S. Food and Drug Administration. The methods and performance characteristics have been reviewed and approved by the New York State Department of Health.

<p>Revision: 1 © This document is copyright of Eurofins</p>	<p>Effective date: May 18, 2018 During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 2 of 8</p>
--------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------	----------------------------------

Test Requisition Instructions



Prenatal Screening Test Requisition Form Instructions

- 1 Account Information – Please enter Ordering Physician name and Referring Ob/Gyn name and phone number, if applicable. A provider signature is required for patients with Medicaid.
- 2 Specimen Labels – Preprinted with the requisition number. Please enter the patient's last and first name EXACTLY as they appear on the requisition form. Affix label(s) to patient specimen(s). Please complete date drawn and drawn-by fields.
- 3 Patient Information – For all tests, please complete patient weight, ethnicity and current pregnancy information. Complete additional patient history as appropriate for test(s) ordered. *Please note, all patient information requested is used to ensure the most accurate risk assessment possible for your patient.*
- 4 Gestational Age – Complete for tests other than First Trimester Screen I FB, Sequential Screen I FB, PreeclampsiaScreen™ | T1 or Maternal Fetal Screen™ | T1 which require CRL (see section 6).
- 5 Biophysical Information – Complete this section for preeclampsia screening only.
- 6 Ultrasound Information – Please provide sonographer and supervisor names and credentialing numbers. Enter all ultrasound information as appropriate for test(s) ordered.
- 7 Test Requests – Tests are ordered by specimen type. Check all tests that apply and provide appropriate ICD codes.
- 8 Cell Free DNA – BOTH the physician and patient signatures are required.
- 9 Billing Information – Provide photocopy of front and back of insurance card or print the information in the required fields.
- 10 Patient Signature – Required for all tests.

1 Account Information – Please enter Ordering Physician name and Referring Ob/Gyn name and phone number, if applicable. A provider signature is required for patients with Medicaid.

2 Specimen Labels – Preprinted with the requisition number. Please enter the patient's last and first name EXACTLY as they appear on the requisition form. Affix label(s) to patient specimen(s). Please complete date drawn and drawn-by fields.

3 Patient Information – For all tests, please complete patient weight, ethnicity and current pregnancy information. Complete additional patient history as appropriate for test(s) ordered. *Please note, all patient information requested is used to ensure the most accurate risk assessment possible for your patient.*

4 Gestational Age – Complete for tests other than First Trimester Screen I FB, Sequential Screen I FB, PreeclampsiaScreen™ | T1 or Maternal Fetal Screen™ | T1 which require CRL (see section 6).

5 Biophysical Information – Complete this section for preeclampsia screening only.

6 Ultrasound Information – Please provide sonographer and supervisor names and credentialing numbers. Enter all ultrasound information as appropriate for test(s) ordered.

7 Test Requests – Tests are ordered by specimen type. Check all tests that apply and provide appropriate ICD codes.


8 Cell Free DNA – BOTH the physician and patient signatures are required.

9 Billing Information – Provide photocopy of front and back of insurance card or print the information in the required fields.

10 Patient Signature – Required for all tests.

Please call 1-888-NTD-LABS (683-5227) for further assistance.

Copyright ©2017, Eurofins NTD, LLC. All rights reserved. NTD is a registered trademark of Eurofins Scientific. All other trademarks are the property of their respective owners. NTD-51111-0817 Printed in USA

	<p align="center">Document name: GEN-2-043 Attachment N-Amniotic Fluid AChE Only Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059815 NTD Labs SOP ID: GEN-2-043 (Att N) Revision:1</p>
-----------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Collection of Amniotic Fluid Specimens:

Collect Amniotic Fluids specimens following your standard operating procedure. Centrifuge to remove cells and transfer supernatant to a transfer tube. Send the supernatant for AFP testing.

Special handling needs between time of collection and time received by laboratory

Do not freeze specimens. Label all specimen tubes. Close shipping canisters tightly. Ship specimens within 24 hours via FedEx (priority overnight and Saturday delivery). Transport at room temperature. Refrigerate specimens if delayed before shipping.

Positive Patient Identification and Specimen Labeling

The clinician collecting the specimen is responsible for correctly identifying the patient using two unique patient identifiers that include the patient's complete first and last name, date of birth and/or medical record or hospital number.

Every patient must verbalize his/her name to the clinician, if able to do so by asking "Would you please tell me (or spell) your name and birthdate". It is unacceptable for the clinician to ask the patient to confirm his/her name that was verbalized by the clinician.

The Specimen MUST be labeled with two forms of patient ID that match EXACTLY with the information on the test requisition. The test Requisition has peel off labels preprinted with the Test Requisition Number that may act as one form of ID. Please print the patient's full first and last names on the label exactly as they appear on the requisition and affix to the specimen. Please retain the requisition number in your records and use if inquiring about test results.

If using your own specimen labels please include two unique patient identifiers that include the patient's complete first and last name, date of birth and/or medical record or hospital number and/or the test requisition number. All information must match exactly on the test requisition form.

Enter the name of the person drawing the specimen on the test requisition to allow for traceability

Safe disposal of materials used in collection


Use only an approved disposal container for syringes and needles used in collection of specimens. Do not throw needles or syringes into the garbage, recycling containers, the toilet, plastic milk jugs, bleach bottles, or soda bottles. NTD cannot accept biohazardous materials and is not affiliated with any licensed medical waste haulers. If necessary please contact a regional hospital to direct you to an appropriate medical waste removal expert.

Revision: 1	Effective date: May 18, 2018	Page 4 of 8
© This document is copyright of Eurofins		During the actual use of this form, cover page need not be printed.

NTD Prenatal Screening Requisition Form

NTD Eurofins Clinical Diagnostics 80 Roland Rd, Suite 1 • Melville, NY • 11747 ntdlabs.com • Phone: 855-754-5221		Prenatal Screening Requisition	
Physician information			
Ordering Physician	Ordering Physician Signature	Physician Code	
Referring US/Phys	Referring US/Phys Phone		
Patient Information			
Last Name	First Name	Date of Birth	Medical Record #
Address		City	State Zip Phone
Due Date	<input type="checkbox"/> By LMP <input type="checkbox"/> By US	Weight (lbs)	Height (in)
*Sequential Screen (F) is dated based on the first trimester CRL.			
Ethnicity <input type="checkbox"/> African American or Caribbean <input type="checkbox"/> Ashkenazi Jewish <input type="checkbox"/> Asian <input type="checkbox"/> Asian Indian <input type="checkbox"/> Caucasian <input type="checkbox"/> Hispanic <input type="checkbox"/> Native American <input type="checkbox"/> Other			
Current Pregnancy (check all that apply) <input type="checkbox"/> IVF - Age of Egg at Harvest _____ yrs. <input type="checkbox"/> Twin <input type="checkbox"/> Multiple # _____ <input type="checkbox"/> Smoker			
Pregnancy History (check all that apply) <input type="checkbox"/> Prior Pregnancy with Down syndrome <input type="checkbox"/> Prior pregnancy with Trisomy 13 <input type="checkbox"/> Prior Pregnancy with Trisomy 18			
ONTD History (check all that apply) <input type="checkbox"/> Family Hx of ONTD (relationship to patient) _____ <input type="checkbox"/> Valproic Acid (Depakene) or Carbamazepine (Tegretol) THIS Pregnancy <input type="checkbox"/> Insulin dependent Before Pregnancy _____ <input type="checkbox"/> Different address than above during the first 3 months of pregnancy.			
Preeclampsia History (check all that apply) <input type="checkbox"/> Previous Pregnancy with Preeclampsia <input type="checkbox"/> Previous delivery > 24 weeks <input type="checkbox"/> Patient's mother with history of Preeclampsia <input type="checkbox"/> History of Chronic Hypertension			
Biophysical Information (for Preeclampsia Screen T1)			
Height (ft) (in)	Blood Pressure Date	Left Arm Blood Pressure	Right Arm Blood Pressure
Ultrasound Information			
Sonographer	FMB or NIDEP	Sonographer to provider	FMB or NIDEP
Ultrasound Date	CRL (43-84mm) _____ mm	NT _____ mm	NB <input type="checkbox"/> Present <input type="checkbox"/> Absent
Twin B	CRL (43-84mm) _____ mm	NT _____ mm	NB <input type="checkbox"/> Present <input type="checkbox"/> Absent
<input type="checkbox"/> Monochorionic <input type="checkbox"/> Dichorionic		PE Risk Not Calculated in Twins	
First Trimester Test Requests			
<input type="checkbox"/> Maternal Fetal Screen T1 (PIGF, AFP, PAPP-A, free Beta, Inhibin-A, NT w/optional NB) (10w0d - 13w6d)		Serum - SST (Red/Grey or Gold Top) or Red Top	NTD Labs Use Only
<input type="checkbox"/> Preeclampsia Screen T1 (PIGF, AFP, PAPP-A w/ optional UTAD and MAP) (10w0d - 13w6d)		Serum - SST (Red/Grey or Gold Top) or Red Top	NTD Labs Use Only
<input type="checkbox"/> First Trimester Screen Fβ (Free Beta, PAPP-A, AFP, NT w/optional NB) (9w0d - 13w6d)		Dried Blood Spot	NTD Labs Use Only
<input type="checkbox"/> Cystic Fibrosis Carrier Screening <input type="checkbox"/> Male (please provide female reproductive partner) Female Name _____ Female DOB _____		Dried Blood Spot	NTD Labs Use Only
Second Trimester Test Requests			
<input type="checkbox"/> Sequential Screen Fβ (free-Beta, AFP, uE3, Inhibin-A + First trimester Screen) (15w0d - 21w6d) <i>Patient must have first trimester screen done through NTD Labs</i> First Trimester Patient ID Number _____		Serum - SST (Red/Grey or Gold Top) or Red Top	NTD Labs Use Only
<input type="checkbox"/> Quad Screen Fβ (free-Beta, AFP, uE3, Inhibin-A) (15w0d - 21w6d)		Serum - SST (Red/Grey or Gold Top) or Red Top	NTD Labs Use Only
<input type="checkbox"/> AFP Test (for ONTD) (15w0d - 21w6d) <input type="checkbox"/> Repeat test for Elevated MSAFP		Serum - SST (Red/Grey or Gold Top) or Red Top	NTD Labs Use Only
<input type="checkbox"/> Second Trimester Screen Fβ (Free Beta, AFP) (15w0d - 21w6d)		Dried Blood Spot	NTD Labs Use Only
<input type="checkbox"/> AFP Test (for ONTD) (15w0d - 21w6d) <input type="checkbox"/> Repeat test for Elevated MSAFP		Dried Blood Spot	NTD Labs Use Only
Amniotic Fluid Specimen Test Requests			
<input type="checkbox"/> AF-AFP with reflexive AChE (15w0d - 21w6d)		Amniotic Fluid	NTD Labs Use Only
<input type="checkbox"/> Amniotic Fluid AChE Only		Amniotic Fluid	NTD Labs Use Only
Billing Information (Please Attach a Copy of The Front and Back of The Patient's Insurance Card or Provide Information Below)			
Insurance Company	Plan Name	Group #	
Subscriber's Last Name, First Name	Insurance ID#	Referral Authorization #	
Insurance Claims Address	Secondary Insurance Information		
<small>I authorize Eurofins NTD, LLC to obtain and release relevant medical and other information and to directly bill and submit claims to Medicare, Medicaid, Medicare Supplemental and/or insurance providers ("insurance") for laboratory/ medical services that Eurofins NTD, LLC provides to me. I assign insurance benefits to Eurofins NTD, LLC and acknowledge that charges that are not covered by insurance, including any applicable co-payments, deductibles, co-insurance, non-covered charges, and charges due to no authorization are my responsibility and I agree to pay for such charges.</small>			
Patient Signature (Required for all Tests)			Date
Specimen Labeling			
Date Drawn		Drawn By:	
<div style="border: 1px solid black; padding: 5px; display: inline-block;"> UR4306001 Last Name, First Name </div>		<div style="border: 1px solid black; padding: 5px; display: inline-block;"> UR4306001 Last Name, First Name </div>	
<small>Enter patient's name on specimen identification label(s) EXACTLY as it appears on the Requisition Form below Two forms of patient ID MUST appear on both the Test Requisition Form and the specimen</small>			

NTD-51101-0118


	<p align="center">Document name: GEN-2-043 Attachment N-Amniotic Fluid AChE Only Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059815 NTD Labs SOP ID: GEN-2-043 (Att N) Revision:1</p>
-----------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

EMR Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Required
SECTRIACHE	HXNTD	Family History of ONTD <ul style="list-style-type: none"> • None • Prev Child • Patient's Sister's Child • Patient's Brother's Child • Patient's Sister • Patient's Brother • Patient's Aunt • Patient's Uncle • Partner's Prev Child • Partner's Brother 	Answer List	YES
SECTRIACHE	NOF	Number of Fetuses <ul style="list-style-type: none"> • 1 • 2 	Answer List	YES
SECTRIACHE	EDD	Due Date	Plain Text	NO
SECTRIACHE	GAUS	GA at U/S (Weeks.Days)	Plain Text	NO
SECTRIACHE	USDATE	U/S Date	Plain Text	NO

EMR Result Codes:


Data Type	Code	LOINC	Name	Contains Result	Comments
CE	MKRANA		Markers/Analytes	No	Included if any 1T markers are available
CE	2T		2nd Trimester	No	
CE	ACHE	30106-9	Acetylcholinesterase	Yes	Measurements contained in NTEs
CE	DGD		Demographic Data	Yes	Included if demographic data is available (contained in NTEs)
CE	2TD		2 nd Trimester Data	Yes	Included if test-specific data is available (contained in NTEs)
CE	REC		Recommendations	Yes	Only displayed if recommendations are available
CE	COM	55107-7	Comments	Yes	Only displayed if comments are available
CE	FTR		Footer	Yes	Only displayed if footer is available
CE	NOT		Notification	Yes	Included for Unsatisfactory Specimens Only

	<p align="center">Document name: GEN-2-043 Attachment N-Amniotic Fluid AChE Only Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059815 NTD Labs SOP ID: GEN-2-043 (Att N) Revision:1</p>
-----------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Example Report

Ache Only Report

<p>Physician ID #: 24328 Physician Tel #: (000) 000-0000 OB SPECIALISTS 100 ANYWHERE ST MELVILLE, NY 11747</p>	<p>Patient Name: ACHE, ONLY Client ID #: Patient ID #: 17AF1101172 <u>2nd Trimester data</u> Date of Birth: 10/08/88 Insulin Rx: Yes Age at EDC: 29 Fam Hx ONTD: No Mat. Weight: State: NY Ethnicity: Afr. Amer./Carib. Draw Date: 10/31/17 Prev Chrom Hx: None GA @ Draw: 23w2d Multi. Preg: No GA by: EDC By U/S Date Received: 11/01/17 EDC: 02/25/18 Report Date: 11/03/17</p>
------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



2nd Trimester


Marker/Analyte	Value
Acetylcholinesterase	Negative

Open Neural Tube Defect NOT suspected.

Comments:


- It is the experience of this program that amniotic fluids from pregnancies affected with ONTDs produce positive AChE results. Malformations of the abdominal wall of the fetus (gastroschisis or omphalocele) produce positive AChE in approximately 80% of affected cases. AChE results have been negative in cases of congenital nephrosis.

 Jonathan B. Carmichael, Ph.D
Laboratory Director,
Eurofins NTD, LLC

 Terrence W. Hallahan, Ph.D
Laboratory Director,
Eurofins NTD, LLC

Physician Information Brochure

<p>Revision: 1</p>	<p>Effective date: May 18, 2018</p>	<p align="right">Page 7 of 8</p>
<p>© This document is copyright of Eurofins During the actual use of this form, cover page need not be printed.</p>		

	<p align="center">Document name: GEN-2-043 Attachment N-Amniotic Fluid AChE Only Instruction Manual</p>	<p>Eurofins Document Reference: 1-D-QM-CF -9059815 NTD Labs SOP ID: GEN-2-043 (Att N) Revision:1</p>
-----------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------

N/A

Patient Information Brochure

N/A

<p>Revision: 1 © This document is copyright of Eurofins</p>	<p>Effective date: May 18, 2018 During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 8 of 8</p>
------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------	------------------------------------------------