

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	--	--

GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual

Eurofins Document Reference	1-D-QM-CF -9059823	Type of document	CF - Controlled Form
NTD Labs SOP ID	GEN-2-043 (Att E)	Division	1-D Clinical Diagnostics Services
Effective Date	Aug 1, 2018	Business Line(s) / Unit(s)	(1-DU) Clinical Diagnostics US : (2-59) Clinical Diagnostics Services North-East US : (EUUSME2) Eurofins NTD (US)
Status	Effective	Periodic Review Date	Jul 31, 2020
		Functional Area	QM - Quality Management

Written by	Elizabeth Sylander
Functional Document Owner	Eurofins D CDS US Laboratory Management; Eurofins D CDS US Quality Management Department; Eurofins D CDS US Quality Control Department
Review and Approval	<ul style="list-style-type: none"> Reviewers: Stephanie Zichi; Margaret Palladino; David Krantz; Norman Moore; Christina Deer; Lisa Schmitt; Jonathan Hayden Approver (Laboratory Director Only): Terrence Hallahan
Reason for Revision	Remove specimen collection section

Revision Log

Date	Rev.	Author	Description
May 18, 2018	1	Eurofins D CDS US Laboratory Management; Eurofins D CDS US Quality Management Department; Eurofins D CDS US Quality Control Department	

Electronic Signatures

Stephanie Zichi;Review;Jul 3, 2018 9:24 AM EDT Christina Deer;Review;Jul 3, 2018 9:31 AM EDT David Krantz;Review;Jul 5, 2018 9:39 AM EDT Jonathan Hayden;Review;Jul 6, 2018 4:56 PM EDT Margaret Palladino;Review;Jul 9, 2018 11:27 AM EDT Norman Moore;Review;Jul 10, 2018 8:32 AM EDT Lisa Schmitt;Review;Jul 11, 2018 10:30 AM EDT Terrence Hallahan;Approval;Jul 18, 2018 9:52 AM EDT
--

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="center">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	---

Reporting Title: Verifi®

Test Definition: CFDNA1
Testing Location: Redwood City, CA
Reporting Location: Melville, NY

Description:

Circulating cell-free DNA is examined from maternal whole-blood to identify fetal aneuploidies Trisomy 21, Trisomy 18 and Trisomy 13. Optionally, the test may be used to identify sex chromosome abnormalities (monosomy X, XXX, XXY, and XYY), Aneuploidy (trisomy) of all chromosomes, including sex chromosome aneuploidies and/or microdeletions such as 1p36 deletion, 4p-(Wolf-Hirschhorn syndrome), 5p- (cri-du-chat syndrome), 15q11.2 (Prader-Willi syndrome/ Angelman syndrome), and 22q11.2 deletion (DiGeorge).

Analytical Method(s):

Nucleic Acid extraction, DNA sequencing, and analysis of sequencing results to determine fetal aneuploidy and fetal microdeletions (if applicable). Test performed at Illumina, Inc. 800 Saginaw Drive, Redwood City, CA 94063.

Patient preparations:

Counsel patient on prenatal screening for detectable genetic abnormalities with cell-free DNA testing.

Specimen Requirements:

Special Timing: Draw blood after 10 weeks, 0 days

Container/Tube: Streck Tube

Specimen Volume: Fill tube to capacity

Specimen Stability: Whole blood specimens are stable at ambient temperature for 5 days.

Specimen Rejection Criteria: hemolysis, lipemia, incorrect tube type

Specimen Collection Instructions

See Blood Specimen Collection from Venipuncture Instruction Manual.

Additional Information:


1. Indications for Testing: Patient is at high risk for aneuploidy. Indications of high risk include advanced maternal age, positive serum screen, abnormal ultrasound and history suggestive of increased risk for the specified chromosomal abnormalities.

CPT Codes: 81420

Reference Values:

Aneuploidy:
Aneuploidy suspected (borderline value)
Aneuploidy detected
No aneuploidy detected

Revision: 2	Effective date: Aug 1, 2018	Page 1 of 16
© This document is copyright of Eurofins		
During the actual use of this form, cover page need not be printed.		

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	--

Sex Chromosomes:

Aneuploidy detected

No aneuploidy detected - Results consistent with two sex chromosomes (XX)

No aneuploidy detected - Results consistent with two sex chromosomes (XY)

Twin pregnancy:

Detected - Results consistent with the presence of Y chromosome material

Not Detected - Results consistent with the absence of Y chromosome material

Microdeletions:

Abnormality detected

Not detected

No microdeletions detected

An Interpretive Report will be provided.

Supplemental Report:

No

Testing Algorithm:

Follow up testing:


1. Aneuploidy Detected or Suspected: Counsel patient about invasive diagnostic procedures such as CVS or amniocentesis.
2. Microdeletion Detected or Suspected: Counsel patient about invasive diagnostic procedures such as CVS or amniocentesis with microarray analysis.

Consents/Authorizations:

Order Practitioners Acknowledgment

The ordering physician MUST have on file at Eurofins NTD a signed 'Acknowledgement of Practitioner's Obligation to Obtain Informed Consent for Genetic Testing' for a test request to be processed.

<p>Revision: 2</p> <p>© This document is copyright of Eurofins</p>	<p>Effective date: Aug 1, 2018</p> <p>During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 2 of 16</p>
--	---	-----------------------------------


	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	--

Acknowledgement of Practitioner's Obligation to Obtain Informed Consent for Genetic Testing

- I understand that The New York State Department of Health, Clinical Laboratory Standards of Practice requires NY-permitted laboratories to notify ordering practitioners that they are required to obtain informed consent for genetic testing from their patients.
- I am aware that Eurofins, NTD, LLC (PFI 3173) ("Eurofins") provides a genetic test consent form that includes test-specific information to aid me in fulfilling my obligations to obtain informed consent.
- I understand that informed consent for the NIPT tests must be obtained utilizing the form specific for that test, copies of which have been provided to me by Eurofins
 - I have been apprised that I may obtain these forms in hard copy from Eurofins, and that the forms are also available on the Eurofins website (http://www.ntdlabs.com/providers/consent_forms.php).
- I acknowledge that this information has been made available to me for patient use in decision-making and the informed consent process.
- I acknowledge that I have reviewed and understand the informed consent for genetic testing form relevant to the testing I am ordering as described above and that the information contained in the form shall be conveyed to each patient and/or their guardian in obtaining full and effective informed consent for each genetic test I order.
- I confirm that I shall maintain documentation that I obtained informed consent for such testing in each of my patients' medical charts.
- I confirm that I will make available and provide a copy of the signed patient consent form to Eurofins NTD, LLC upon request.

Practitioner Signature	Date
Printed Name	
Hospital, Facility or Clinic Name	
Address	
City, State, Zip	
Phone #	
Fax#	
NTD Account #	

Please Return to:
Eurofins NTD, LLC, 80 Ruland Road, Suite 1, Melville, NY 11780
Fax: 631-425-0864

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	--

Patient Informed Consent

The patient **MUST** sign a verifi 'Patient Informed Consent' form that is maintained in the ordering physicians office.



Patient informed consent

INTRODUCTION: This form describes the benefits, risks, and limitations of this screening test. You should seek pre-test counseling by a genetic counselor or other experienced health care provider prior to undergoing this test. Read this form carefully – and ask any questions you may have of your health care provider -- before making your decision about testing.

PURPOSE: The purpose of the Verifi™ Test and the Verifi™ Plus Test is to screen your pregnancy for certain chromosomal abnormalities, also known as "aneuploidies." Both tests give information about whether there may be extra copies (trisomy) of chromosomes 21, 18, and 13, and the option to know if there is an extra copy of a sex chromosome (X or Y), and/or a missing copy of sex chromosome (thMX). Fetal sex may also be reported. The Verifi™ Plus Test has the option to screen for aneuploidies (extra copies) in all chromosomes. In addition, the option to screen for the following microdeletions (small, missing parts of chromosomes) syndromes: 1p36 deletion, 4p- (Wolf-Hirschhorn syndrome), 5p- (cri-du-chat syndrome), 15q11.2 (Prader-Willi syndrome/Angelman syndrome), 22q11.2 deletion (DiGeorge syndrome or velocardiofacial syndrome) is also available. For chromosomes 21, 18, and 13, the Verifi™ Test is validated in singleton and twin pregnancies. In twin pregnancies, sex chromosome testing can only screen for the presence or absence of the Y chromosome, and not for extra or missing sex chromosomes. Both Verifi™ and Verifi™ Plus can be performed as early as 10 weeks 0 days gestational age. Consult your health care provider if you would like more information about this screening test, including risks, limitations, performance data, error rates, descriptions of the conditions being screened, and what these results may mean to your pregnancy.


HOW THIS TEST WORKS: This test screens for specific chromosomal abnormalities by looking at the DNA (genetic material) in your blood. The sample of blood includes a combination of both your DNA and the DNA from the pregnancy. A technology called massively parallel sequencing is used to count the amount of DNA from each test chromosome and/or from specific regions of chromosomes. The laboratory then uses an analysis method to determine if each of the conditions you have elected to test for is likely to be present or absent.

SEX OF PREGNANCY: Depending upon the option you and your health care provider elect, the test results may include the sex of the pregnancy. If you do not wish to know the sex, please tell your health care provider not to disclose this information to you. Depending upon the test ordered, you may not be able to prevent learning the sex of your pregnancy. In rare instances, incorrect sex results can occur.

LIMITATIONS OF THE TEST: These are screening tests that look only for specific chromosomal abnormalities. This means that other chromosomal abnormalities may be present and could affect your pregnancy. A "No Aneuploidy Detected" result does not guarantee a healthy pregnancy or baby and does not eliminate the possibility that your pregnancy may have birth defects, genetic conditions, or other conditions, such as open neural tube defects or autism.

There is a small possibility that the test results might not reflect the chromosomes of the fetus, but may reflect chromosomal changes of the placenta (confined placental mosaicism, CPM) or of you (maternal chromosomal abnormalities). While these tests are not designed to assess your health, in some cases, information about your health may be revealed directly or indirectly (e.g., when combined with other information). Examples include maternal XXX, sex chromosome status or benign or malignant maternal neoplasms. In a twin pregnancy, the status of each individual fetus cannot be determined.

<p>Revision: 2 © This document is copyright of Eurofins</p>	<p>Effective date: Aug 1, 2018 During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 4 of 16</p>
---	--	-----------------------------------

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	--	--



These tests, like many tests, have limitations, including false negative and false positive results. This means that the chromosomal abnormality being tested for may be present even if you receive a negative result (this is called a 'false negative'), or that you may receive a positive result for the chromosomal abnormality being tested for, even though the abnormality is not actually present (this is called a 'false positive').

In the case of a twin pregnancy, the presence or absence of Y chromosome material can be reported. The occurrence of sex chromosome aneuploidies cannot be evaluated in twin pregnancies. In the case of a vanishing twin, the test result may reflect the DNA of the vanishing twin, leading to a higher probability of false positive or false negative results.

No irreversible clinical decisions should be made based on these screening results alone. If definitive diagnosis is desired, chorionic villus sampling or amniocentesis would be necessary. In some cases, other testing may also be necessary. Some rare chromosomal aneuploidies may only occur in mosaic form. Clinical consequences depend on the chromosome involved and cannot be predicted prenatally.

Consult your health care provider for more information about your results and what they may mean for your pregnancy, what options you will have for further testing, and whether additional testing is recommended for you based on your clinical history.

TEST PROCEDURE: A tube of your blood will be drawn and sent to Verinata Health, Inc., a wholly owned subsidiary of Illumina, Inc., which will then analyze your blood.

PHYSICAL RISKS: Side effects of having blood drawn are uncommon, but may include dizziness, fainting, soreness, bleeding, bruising, and, rarely, infection.


DISCRIMINATION RISKS: Genetic information could be used as a basis of discrimination. To address concerns regarding possible health insurance and employment discrimination, some countries, U.S. states, and the U.S. government have enacted laws to prohibit genetic discrimination in those circumstances. The laws may not protect against genetic discrimination in other circumstances, such as when applying for life insurance or long-term disability insurance. Talk to your health care provider or genetic counselor if you have concerns about genetic discrimination prior to testing.

PREGNANCY OUTCOME INFORMATION: Collecting information on your pregnancy after testing is part of a laboratory's standard practice for quality purposes and is required in several states. As such, Illumina or its designee may contact your health care provider to obtain this information. By executing this informed consent, you agree to allow your health care provider to provide this information to Illumina or its designee.

SECONDARY FINDINGS: In the course of performing the analysis for the indicated tests, information regarding other chromosomal alterations, also known as "secondary findings" may become evident. Our policy is to NOT REPORT on any secondary findings that may be noted in the course of analyzing the test data.

PRIVACY: Test results are kept confidential. Your test results will only be released in connection with the testing service, to your health care provider, his or her designee, other health care providers involved in your medical care, or to another health care provider as directed by you (or a person legally authorized to act on your behalf) in writing, or otherwise as required or authorized by applicable law.

<p>Revision: 2 © This document is copyright of Eurofins</p>	<p>Effective date: Aug 1, 2018 During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 5 of 16</p>
--	---	-----------------------------------

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	--



CROSS-BORDER DATA TRANSFER: If you are from outside the United States, your specimen and associated health information will be sent to the United States in order for the testing to be completed. As part of the testing, additional health information about you will be created and maintained. Your country may consider the legal privacy protections in the United States to be inadequate.

USE OF INFORMATION AND LEFTOVER SPECIMENS FOR RESEARCH: Eurofins NTD, LLC performs research to help develop and provide safe and effective screening tests and to contribute to advancing biomedical knowledge. Your sample will be discarded within sixty days after taken, unless you specifically give consent. Permission to allow the use of your de-identified sample in research and development studies is entirely voluntary.


TEST RESULTS: Your test results will be sent to the health care provider.

PATIENT CONSENT: By signing this form, I, the patient having this screening performed, acknowledge that: (i) I have been offered the opportunity to ask questions and discuss with my health care provider the benefits, risks, and limitations of the test to be performed; (ii) I have discussed the test limitations (reliability of positive and negative test results; the predictive value of the test results; and that the test is not a diagnostic test, but a screening test and is not definitive) with the health care provider who ordered the test; (iii) I have been informed about the availability and importance of genetic counseling and have been provided with information identifying an appropriate health care provider from whom I might obtain such counseling; (iv) I have received, read, and understood the Patient Informed Consent in its entirety and that I may retain a copy for my records (v) I consent to having this test performed; and (vi) I will discuss the results and appropriate medical management with my health care provider.

Patient Signature:

Date:


<p>Revision: 2 © This document is copyright of Eurofins</p>	<p>Effective date: Aug 1, 2018 During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 6 of 16</p>
---	--	-----------------------------------

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="center">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	---

Disclaimer:

This test was developed by, and its performance characteristics were determined by Verinata Health, Inc., a wholly-owned subsidiary of Illumina, Inc., a CLIA-certified and CAP-accredited clinical testing laboratory. It has not been cleared or approved by the U.S. Food and Drug Administration.

<p>Revision: 2 © This document is copyright of Eurofins</p>	<p>Effective date: Aug 1, 2018 During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 7 of 16</p>
--	--	-----------------------------------

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="center">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	---

Test Requisition Instructions

Complete Physician Information, Patient Information and Clinical Information sections. Be sure to check off the specific Test Request.

Specimen Labels – Preprinted with the requisition number. Please enter the patients last and first names EXACTLY as they appear on the requisition form. Affix the label to the patient specimen. Please complete the date drawn and drawn-by fields.

Billing Information – Provide a photocopy of the front and back of insurance card or print the information in the required fields. Please provide credit card information to cover tests ordered requiring additional charges.

Special handling needs between time of collection and time received by laboratory

Do not freeze specimens. Label all specimen tubes. Close shipping canisters tightly. Ship specimens within 24 hours via FedEx (priority overnight and Saturday delivery). Transport at room temperature. Refrigerate specimens if delayed before shipping.


<p>Revision: 2 © This document is copyright of Eurofins</p>	<p>Effective date: Aug 1, 2018 During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 8 of 16</p>
---	--	-----------------------------------

Test Requisition

Verifi™
Prenatal Test

Physician Information			
Ordering Physician	Ordering Physician Signature	Physician Code	
Referring OB/Gyn	Referring OB/Gyn Phone		
Patient Information			
Last Name	First Name	Date of Birth ____/____/____	Medical Record #
Address		City	State Zip
Phone		email	
Clinical Information			
Gestational age: Weeks: _____ Days: _____	As estimated on: ____/____/____	DATING METHOD (Must choose one): <input type="checkbox"/> LMP <input type="checkbox"/> Date of implantation <input type="checkbox"/> CRL	
Maternal height: <input type="checkbox"/> cm or <input type="checkbox"/> ft in	Maternal Weight: <input type="checkbox"/> lbs or <input type="checkbox"/> kg	Age of egg donor: _____	
TEST INDICATIONS (Choose at least one): <input type="checkbox"/> Advanced maternal age (≥35 years) <input type="checkbox"/> Positive serum screen <input type="checkbox"/> Abnormal ultrasound <input type="checkbox"/> High suggestive of increased risk for the specified chromosome aneuploidies <input type="checkbox"/> Low risk/maternal anxiety <input type="checkbox"/> Other		Comments: _____	
Test Request			
<input type="checkbox"/> Verifi™ Prenatal Test for Singleton Pregnancy for chromosomes 21, 18, 13 (≥10w0d) <input type="checkbox"/> Sex chromosome aneuploidies			
Verifi™ Plus Prenatal Test for Singleton Pregnancy for chromosomes 21, 18, 13 and additional tests as selected below (≥10w0d) <input type="checkbox"/> Expanded autosomal trisomies* (all chromosomes including sex chromosome aneuploidies) <input type="checkbox"/> Microdeletion panel† [22q11 (DiGeorge), 15q11 (Angelman/Prader-Willi), 1p36, 4p-(Wolf-Hirschhorn), 5p-(Cri-du-chat)] <input type="checkbox"/> Sex chromosome aneuploidies *Additional cost for this test option, please complete the credit card info. section below.			
<input type="checkbox"/> Verifi™ Prenatal Test for Twin Pregnancy for chromosomes 21, 18, 13(≥10w0d) <input type="checkbox"/> Presence of Y chromosome			
Patient Consent for Research			
Eurofins NTD, LLC performs research to help develop and provide safe and effective screening tests and to contribute to advancing biomedical knowledge. My sample will be discarded within sixty days after taken, unless I specifically consent below. Permission to allow the use of my de-identified sample in research and development studies is entirely voluntary.			
<input type="checkbox"/> I give permission for Eurofins NTD, LLC to retain any remaining de-identified sample for future research and development. <input type="checkbox"/> I do not give permission to use the remaining sample in any research studies.			
Patient Signature _____		Date ____/____/____	
Billing Information			
ICD Codes: _____			
Please Attach a Copy of The Front and Back of the Patient's Insurance Card or Provide Information Below			
Insurance Company	Plan Name	Group #	
Subscriber's Last Name, First Name	Insurance ID#	Referral Authorization #	
Insurance Claims Address	Secondary Insurance Information		
I authorize Eurofins NTD, LLC to obtain and release relevant medical and other information and to directly bill and submit claims to Medicare, Medicaid, Medicare Supplemental and/or insurance providers ("insurance") for laboratory medical services that Eurofins NTD, Inc. provides. I assign insurance benefits to NTD Labs and acknowledge that charges that are not covered by insurance, including any applicable co-payments, deductibles, co-insurance, non-covered charges, and charges due to no authorization are my responsibility and I agree to pay for such charges.			
Patient Signature _____		Date ____/____/____	
Credit Card Information			
Credit Card No. _____	Expiration ____/____	CW	Authorized Amount: \$ _____
<input type="checkbox"/> MasterCard <input type="checkbox"/> Visa <input type="checkbox"/> American Express <input type="checkbox"/> Other: _____		Patient Signature (required) _____	
Signature of Credit Card Holder (if different than patient) _____		Date ____/____/____	
Specimen Labeling			
Sample Type: Streck Cell Free DNA BCT (Beige/Black Top Tube)		Date Drawn ____/____/____	Drawn By _____
<div style="border: 1px solid black; width: 200px; height: 30px; margin: 0 auto;"></div>		<div style="border: 1px solid black; width: 200px; height: 30px; margin: 0 auto;"></div>	
<div style="border: 1px solid black; width: 80px; height: 30px; margin: 0 auto; text-align: center;">NTD Labs Use Only</div>			
Enter patient's name on specimen identification label EXACTLY as it appears on the Requisition Form below. Peel off label and affix to the blood tube.			

NTD-41131-2227

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	--


EMR Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Required
CFDNA1	TESTXY	Test For Sex Chromosome Anomalies <ul style="list-style-type: none"> No Yes 	Answer List	NO
CFDNA1	TESTOTH	Test For Other Chromosomes <ul style="list-style-type: none"> No Yes 	Answer List	NO
CFDNA1	TESTMCD	Test for Microdeletions <ul style="list-style-type: none"> No Yes 	Answer List	NO
CFDNA1	TESTTWINY	If Twins, Test for Presence of Y Chromosome <ul style="list-style-type: none"> No Yes 	Answer List	NO
CFDNA1	TESTIND1	Test Indication <ul style="list-style-type: none"> Advanced Maternal Age Positive Serum Screen Abnormal Ultrasound Hx Suggestive of increased risk for the specified chromosome abnormalities Low Risk Maternal Anxiety Other 	Answer List	YES
CFDNA1	TESTINDOTH	Other Test Indications (Describe)	Plain text	YES
CFDNA1	GA	Gestational Age Weeks/Days	Plain text	YES
CFDNA1	GADATE	Date of Gestational Age Assessment	Plain text	YES
CFDNA1	MATWT	Maternal Weight (lbs)	Plain Text	NO
CFDNA1	MWLBSKGS	Units <ul style="list-style-type: none"> LBS KGS 	Answer List	YES
CFDNA1	HGHT	Maternal Height (Ft.In)	Plain text	NO


EMR Result Codes:

Data Type	Code	LOINC	Name	Contains Result	Comments
CE	SUM	77011-5	Summary Result	Yes	
CE	CHR21	77012-3	Chromosome 21	Yes	Interpretation in NTE
CE	CHR13	77013-1	Chromosome 18	Yes	Interpretation in NTE
CE	CHR13	75977-9	Chromosome 13	Yes	Interpretation in NTE
CE	SEX	77011-5	Sex Chromosome	Yes	Included, if available, Interpretation in NTE
CE	CHRY	77011-5	Y Chromosome	Yes	Included, if available, Interpretation in NTE, For Twins Only
CE	CHROTH		All Other Autosomal	Yes	Included, if available, Interpretation in

Revision: 2 © This document is copyright of Eurofins	Effective date: Aug 1, 2018 During the actual use of this form, cover page need not be printed.	Page 10 of 16
---	--	---------------

	Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual	Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2
---	---	---

			Chromosomes		NTE
CE	MCD		Microdeletion	Yes	Interpretation in NTE
NM	FFRAC		Fetal Fraction	Yes	Included, if available
CE	COM-C		Clinical Comments	Yes	Only displayed if comments are available
CE	COM-O		Other Comments	Yes	Only displayed if comments are available
CE	DGD		Demographic Data		Included if demographic data is available (contained in NTEs)
CE	NIPTMETH		Test Method		Test Method Listed in NTEs
CE	PERF		Prenatal Test Performance Data		Performance Data listed in NTEs
CE	NIPTLIM		Limitations of Test		Limitations of Test Listed in NTEs
CE	FTR		FOOTER		Only displayed if footer is available
CE	NOT		NOTIFICATION		Included for Unsatisfactory Specimens Only

	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p align="right">Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	--	--

Example Report



80 Roland Rd, Suite 1 - Melville, NY - 11747 • (855) 754-8221 • Fax (631) 425-0854 • reporting@ntd-eurofins.com

verifi® from NTD Labs

Physician ID #: 9999 Physician Tel #: (999) 999-9999 OB/GYN SPECIALISTS SUIT 100 SOMEWHERE, US, 99999	Patient Name: DOE, JANE Client ID #: 1234567 Patient ID #: 16SE5999999 Date of Birth: 09/07/71 Draw Date: 09/07/15 Date Received: 09/09/15 Report Date: 09/22/15
--	---

verifi® from NTD Labs was performed at:

Illumina, Inc. 800 Saginaw Drive, Redwood City, CA 94063

See attached report.

The information contained in this communication is intended only for the user or the individual or entity to whom or to which it is addressed, and may contain information that is privileged, confidential and exempt from disclosure under applicable law. If the reader of the communication is not the intended recipient, or the employee or agent responsible for delivering this communication to the intended recipient, you are hereby notified that any dissemination, distribution, copying or disclosure of this communication is strictly prohibited. If you have received this communication in error, please notify the sender immediately by telephone and fax and destroy any copies. If this communication has been received electronically, delete the original communication and any attachment from any computer, server or other electronic recording or storage device or medium. No confidentiality or privilege is waived or lost by any mistransmission.

Please notify NTD if you are receiving patient information or other confidential information and your fax machine is not in a secured location. Please notify us if there is an alternative fax number we may use.



Jonathan B. Carmichael, Ph.D
Laboratory Director,
Eurofins NTD, LLC



Terrence W. Halahan, Ph.D
Laboratory Director,
Eurofins NTD, LLC



FAX: 9999999999

Page 1 of 1

This report was generated on: 09/22/2015 03:24:23 PM

Revision: 2 © This document is copyright of Eurofins	Effective date: Aug 1, 2018 During the actual use of this form, cover page need not be printed.	Page 12 of 16
---	--	---------------

Verifi Plus™ prenatal screening report

Verifi Plus™ Prenatal Test

REPORT RELEASED
 Date: 09/22/15 Time: 08:45 AM

PROVIDER INFORMATION	SECOND RECIPIENT	PATIENT INFORMATION	SAMPLE INFORMATION
Eurofins NTD, Inc Attn: Terrence Halahan, Ph.D. 80 Buland Road Melville, NY 11747 Phone: (888) 754-5221 Fax: (631) 425-0864		Name: Jane Doe DOB: 09/07/1971 GA: XXXXXXXX Indication: XXXXXXXX Medical record/patient ID: XXXXXXXX	Client Sample ID: 1234567 Order ID: 1234567 Date of Draw: 09/07/15 Date Received: 09/22/15 Pregnancy Type: Singleton

ANEUPLOIDY DETECTED

RESULTS SUMMARY:

CHROMOSOME	RESULTS	PPV (%)
Chromosome 21	POSITIVE: Aneuploidy detected Results consistent with pregnancy at increased risk for trisomy 21.	XX
Chromosome 18	NEGATIVE: No aneuploidy detected Results consistent with two copies of chromosome 18.	
Chromosome 13	NEGATIVE: No aneuploidy detected Results consistent with two copies of chromosome 13.	
All Other Autosomal Chromosomes	NEGATIVE: No aneuploidy detected Results consistent with two copies of all other autosomes.	NA*
Sex Chromosomes	POSITIVE: Aneuploidy detected Results consistent with pregnancy at increased risk for XO (monosomy X), XXX, XXY, XYY.	NA**
Microdeletions (1p36, 4p16.3, 5p15.2, 15q11.2, 22q11.2)	POSITIVE: Abnormality detected Results consistent with a microdeletion in the 22q11.2 region. NEGATIVE: No aneuploidy detected Results consistent with no microdeletions detected in the regions of 1p36, 4p16.3, 5p15.2, 15q11.2, 22q11.2	NA***

CLINICAL COMMENTS: This is a screening test; therefore, false positive and false negative results can occur. Results may be reflective of fetal, placental, or maternal conditions. No irreversible clinical decisions should be made based on these screening results alone. Clinical correlation is indicated. If definitive diagnosis is desired, chorionic villus sampling or amniocentesis would be necessary. Genetic counseling is recommended. The fetal fraction (FF) is estimated to be 23%. FF estimation is one component of the Verifi Plus™ algorithm and is combined with other quality metrics to determine the confidence in the results. The FF estimate is not used in isolation to exclude samples.

Positive predictive value (PPV) is calculated based on stated performance, maternal and gestational age as provided on the Test Requisition Form (TRF). Other factors may impact the patient specific PPV. For more information about PPV please visit us at: www.illumina.com/ppv.

*Aneuploidies involving chromosomes other than 21, 13, 18, X, or Y are rare, and prevalence is not available, which precludes accurate calculation of PPV.

**Performance data for sex chromosome aneuploidy is limited, precluding accurate calculation of PPV.

***Performance data for microdeletion analysis is limited, precluding accurate calculation of PPV.

OTHER COMMENTS: Verifi Plus™ comments will be in this section. *Loeam ipsum dolor sit amet, ad diam repudera ius, sed ridere oportet definitibus, odio munere possit an vel. Vellit legere vim uia, vix maecenas reprehenderit et, sicut blandi detracto eos et. Eit legere vim eu, vix maecenas reprehenderit et, sicut blandi detracto eos et dolor consetetur dispartit uia, an meli aliquam informidans.*

PERFORMANCE AND LIMITATIONS

LIMITATIONS OF THE TEST: The Verifi Plus™ Prenatal Test is validated for aneuploidy of any chromosome, including 21, 13, 18, X, and Y and for specific deletions in chromosomal regions: 1p36, 4p16.3, 5p15.2, 15q11.2, 22q11.2, in singleton pregnancies, with gestational age of at least 10 weeks 0 days. This is a screening test that looks only for specific chromosomal abnormalities. A normal result does not eliminate the possibility that the pregnancy is associated with other chromosomal or subchromosomal abnormalities, birth defects, genetic conditions, or other conditions, such as open neural tube defects or autism.

There is a small possibility that the test results might not reflect the chromosomes of the fetus, but may reflect chromosomal changes of the placenta (confined placental mosaicism, CPM) or of you (maternal chromosomal abnormalities). Examples include maternal XX,X sex chromosome status, or benign and malignant maternal neoplasm. CPM may be associated with a higher chance for pregnancy complications or for uniparental disomy (UPD), which may affect the growth and development of the fetus. Some of these rare chromosomal aneuploidies may only occur in mosaic form. Clinical consequences depend on the chromosome involved and can not be predicted precisely. This test, like many tests, have limitations, including false negative and false positive results. A negative test result does not eliminate the possibility of chromosomal abnormalities for the tested chromosomes or microdeletions. See performance metrics for test performance.

PERFORMANCE METRICS†

Chromosome	N	Sensitivity	95% CI	Specificity	95% CI	Accuracy	95% CI
21	510	99.9% (97.9)	99.0 - 100.0	99.9% (99.9)	99.7 - 100.0	-	-
18	501	97.6% (97.9)	95.2 - 99.9	99.9% (99.9)	99.5 - 100.0	-	-
13	501	97.5% (97.9)	95.1 - 98.5	99.9% (99.9)	99.2 - 100.0	-	-
Chromosome	N	Sensitivity	95% CI	Specificity	95% CI	Accuracy	95% CI
Monosomy X	508	95% (95.0)	75.4-99.9	99.0% (99.9)	97.8-99.7	-	-
XX	508	97.6% (97.6)	96.8-99.1	99.0% (99.0)	97.2-99.9	98.0%	95.9-99.3
YY	508	99.1% (99.1)	98.4-99.9	99.9% (99.9)	99.9-99.9	99.0%	97.7-99.7
XXX/YYY/YY		Other sex aneuploidies will be reported if detected. (Limited data of these more rare aneuploidies preclude performance calculations.)					
Microdeletions & other autosomal aneuploidies		Microdeletions and other autosomal aneuploidies (if requested and detected) will be reported. (Limited data of these more rare abnormalities preclude performance calculations.)					

†Data on file at Illumina, Inc. regarding Performance and Method Comparison studies

TEST METHOD: Nucleic Acid extraction, DNA sequencing, and analysis of sequencing results to determine fetal aneuploidy.


DISCLAIMER: The manner in which this information is used to guide patient care is the responsibility of the health care provider, including advising for the need for genetic counseling or diagnostic testing. Any test should be interpreted in the context of all available clinical findings.

DISCLOSURE: The Verifi Plus test was developed by, and its performance characteristics were determined by Verinata Health, Inc. a wholly owned subsidiary of Illumina, Inc. The VHI laboratory is CAP-accredited and certified under the Clinical Laboratory Improvement Amendments (CLIA) as qualified to perform high complexity clinical laboratory testing. The verifi test has not been cleared or approved by the U.S. Food and Drug Administration.

This prenatal test is performed by Verinata Health, Inc., a wholly owned subsidiary of Illumina, Inc. Illumina, Inc., 5000 Greatwood Drive, Redwood City, CA 94063 | 1-855-205-9550
 CAP: 7516112 CLIA: 0602013951 California License #: CLF00040177

Laboratory Director: Meredith Hulse Miller, MD, William K. Seibow, PhD, FACMG, Ethan De Foa, PhD, Beth Sanderson
 © 2017 Illumina, Inc. All rights reserved. Illumina, Verinata Health, the pumpkin orange color, Verifi™ and Verifi Plus™, are trademarks of Illumina, Inc. and/or its affiliates in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.



	<p align="center">Document name: GEN-2-043 Attachment E-verifi® from Eurofins NTD Instruction Manual</p>	<p>Eurofins Document Reference: 1-D-QM-CF -9059823 NTD Labs SOP ID: GEN-2-043 (Att E) Revision:2</p>
---	---	--

Physician Information Brochure

<p>Revision: 2 © This document is copyright of Eurofins</p>	<p>Effective date: Aug 1, 2018 During the actual use of this form, cover page need not be printed.</p>	<p align="right">Page 14 of 16</p>
--	--	--

Patient Information Brochure



What is Verifi™?

Verifi™ is a simple blood test that screens for the most common chromosomal abnormalities that can affect your baby's future health. A sample can be drawn in your doctor's office as early as the 10th week of pregnancy. Verifi™ is available for both singleton and twin pregnancies.*

What does Verifi™ screen for?

Verifi™ looks for too few or too many copies of chromosomes. Missing or extra copies of chromosomes are called "aneuploidies" and are often associated with intellectual or physical disabilities, with different levels of severity.

The most commonly seen aneuploidies in pregnancy include trisomy 21 (Down syndrome), trisomy 18 (Edwards syndrome), and trisomy 13 (Patau syndrome), all of which can be accurately detected with Verifi™.

If you and your healthcare provider choose, Verifi™ can analyze sex chromosomes in order to tell you the gender of your baby and identify sex chromosome abnormalities.

* The Verifi™ test for chromosomes 21, 18, and 13, is available for singleton and twin pregnancies. Sex chromosomes aneuploidy testing is available for singleton pregnancies only. An optional test for twin pregnancies tests for the presence of the Y chromosome.

About Eurofins NTD

For more than 30 years, Eurofins NTD has pioneered the research and development of prenatal screening protocols for Open Neural Tube Defects, Down syndrome, trisomy 18 and trisomy 13. Today, Eurofins NTD serves universities, medical centers, hospitals, laboratories, obstetricians and maternal fetal medicine specialists worldwide—providing risk assessment services that help healthcare professionals and patients make more informed decisions.

To learn more about Verifi™, please speak with your healthcare provider, or call us at 1-866-NTD-LABS (683-5227).

References

1. Bianchi DW, Platt RW, Goldberg JD, et al. Genome-wide fetal aneuploidy detection by maternal plasma DNA sequencing. *Obstet Gynecol*. 2012;119:890–901.
2. Futch T, Spinosa J, Bhatt S, de Foa E, Iyva RP, Schrier AJ. Initial clinical laboratory experience in noninvasive prenatal testing for fetal aneuploidy from maternal plasma DNA samples. *Prenat Diagn*. 2017;37:568–74.



90 Ruland Road, Suite 1, Melville, NY 11747
 www.ntd-eurofins.com | 1-866-NTD-LABS (683-5227)

THIS TEST WAS DEVELOPED BY, AND ITS PERFORMANCE CHARACTERISTICS DETERMINED BY, ILLUMINA. IT IS NOT CLEARLY APPROVED BY THE U.S. FOOD AND DRUG ADMINISTRATION. Copyright ©2018, Eurofins NTD, LLC. All rights reserved. 18-2018-0418 Printed in USA





How quickly will I get my results?

Your results will be sent to your healthcare provider, who will usually receive them within one week.

What do my test results mean?

Your results will tell your healthcare provider whether or not trisomies 21, 18, 13, or sex chromosome abnormalities (if ordered) are likely to be present in your pregnancy. In the case of a positive result, your healthcare provider will discuss what the results mean to your pregnancy, as well as further testing options to consider.

If the sex chromosome option is ordered, results will be reported as either No Aneuploidy Detected or Aneuploidy Detected.

How does the test work?

A sample of your blood is drawn and the genetic material (DNA) from you and your baby is tested.

Verifi™ takes a deeper approach to the science, using an advanced technology called "Whole Genome Sequencing" to analyze millions of DNA fragments per sample and accurately count the number of chromosomes present.

Is Verifi™ right for me?

This screening test is usually offered to pregnant women who have a confirmed singleton or twin pregnancy of at least 10 weeks gestational age, and meet any of the following criteria:

- You are considered to be of advanced maternal age at time of delivery (35 years or older for a singleton pregnancy or 32 years or older for a twin pregnancy)
- You have an abnormal or "positive" serum screen
- Your ultrasound shows concerns or abnormalities with your baby's growth and/or development
- You have a personal or family history suggestive of trisomies 21, 18, 13, or other sex chromosome abnormalities



What is Verifi™ Plus?

The Verifi™ Plus for singleton pregnancies tests for chromosomes 21, 18, 13 and additional tests that you and your healthcare provider may choose to select, as shown below:

- Sex chromosome aneuploidies (monosomy X, XXX, XXY, and XYY) are included if requested; fetal sex (XX or XY) will be reported if no sex chromosome aneuploidy is detected.
- Aneuploidy (trisomy) of all chromosomes, including sex chromosome aneuploidies.
- Select microdeletions such as 1p36 deletion, 4p- (Wolf-Hirschhorn syndrome), 5p- (cri-du-chat syndrome), 15q11.2 (Prader-Willi syndrome/ Angelman syndrome), and 22q11.2 deletion (DiGeorge).

Do normal Verifi™ test results mean that my baby will be perfectly healthy?

Verifi™ is a highly accurate non-invasive screening test. No test, however, can guarantee a baby will not have any medical issues.

Verifi™ addresses only aneuploidies of chromosomes 21, 18, 13, and sex chromosomes* (if ordered). It does not test for or report all genetic and non-genetic problems that may be present in a baby. If the test result is No Aneuploidy Detected, indicating a negative result, it does not completely rule out all potential problems with chromosomes 21, 18, and 13, or all sex chromosome aneuploidies in your baby.

Genetic counseling before and after testing is recommended. If result is 'Aneuploidy Detected' you're considered positive. Women who receive this result should be offered invasive prenatal procedures for confirmation. A negative test does not ensure an unaffected pregnancy.

* Sex chromosome aneuploidy testing is available for singleton pregnancies only.