

MM530 Cystic Fibrosis: *CFTR* Expanded Mutation Panel

Turn around time: 14 business days from sample receipt at the laboratory.

CPT Codes: 81220

Test Description:

Expanded panel that screens for 142 mutations. Use of next generation sequencing and Sanger sequencing for potentially causative variants and areas of low coverage. Accounts for mutations more common in non-Caucasian ethnic groups, as well as rarer mutations across all ethnic groups.

Indications: General population screening of all preconception or pregnant women. Male partner of woman identified as a carrier.

Testing Timeframe: Any time during pregnancy or preconception.

Testing requirements:

Container/Tube: Lavender-top Vacutainer® tube (EDTA).

Specimen Volume: 5 ml of unspun whole blood.

Specimen Stability: Samples are stable at ambient temperature for 6 days

Shipping Information:

Ship samples to the following address:

Eurofins NTD, LLC.
80 Ruland Road, Suite 1
Melville, NY 11747

Questions?

For more information about NTD prenatal screening, ordering test supplies, or billing support, please call us at 1-888-NTD-LABS (683-5227) or email us ntdclientservices@eurofins.com