

MM580 ACOG/ACMG Carrier Screen: Targeted Mutation Panel

Turn around time: 14 business days from sample receipt at the laboratory.

CPT Codes: 81200, 81209, 81220, 81242, 81243, 81251, 81255, 81260, 81290, 81330, 81400

Test Description:

Developed to support both ACOG and ACMG* guidelines and recommendations, the ACOG/ACMG Carrier Screen tests for Cystic Fibrosis, Spinal Muscular Atrophy and 8 other disorders common to those of Ashkenazi Jewish descent.

Indications: General population screening of all preconception or pregnant women. Male partner of woman identified as a carrier of any recessive disorder.

Testing Timeframe: Any time during pregnancy or preconception.

Testing requirements:

Container/Tube: Lavender-top Vacutainer® tube (EDTA).

Specimen Volume: 5 ml of unspun whole blood.

Specimen Stability: Samples are stable at ambient temperature for 6 days

Shipping Information:

Ship samples to the following address:

Eurofins NTD, LLC.
80 Ruland Road, Suite 1
Melville, NY 11747

Questions?

For more information about NTD prenatal screening, ordering test supplies, or billing support, please call us at 1-888-NTD-LABS (683-5227) or email us ntdclientservices@eurofins.com